



CENTRIFUGES

EBA 270

BLOOD / URINALYSIS / PEDIATRIC

Simple. Compact. Complete.

Built-in, horizontal rotor for flat separation of all standard blood, gel, urine and pediatric tubes. Ideal for physicians' and veterinary offices.

MAX. CAPACITY	6 x 15 mL
MAX. RPM/RCF	4,000/2,254

Distributed by:



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.


DIMS. (HxWxD): 9.4" x 12.8" x 15.3"

hettweb.com/eba270

HETTICH EBA 270 APPLICATION PACKAGES



EBA 270 COMPLETE PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO. 270COMPLETE
	6	0.5 Pediatric	10.7 x 36	≤ 4,000	≤ 1,359	
	6	1.6 - 7	13/16 x 75	≤ 4,000	≤ 1,807	
	6	4 - 10	13/16 x 100	≤ 4,000	≤ 2,254	
	6	12 KOVA	17 x 102	≤ 4,000	≤ 2,254	
Includes: (1) Centrifuge, (1) Built-in Horizontal Rotor, (6) Adapters 13s, (6) Adapters 16s, (6) Adapters Pediatric						

Supports over 20 application configurations, please inquire for more information



SMALL IN SIZE, BIG ON RESULTS
 Blood separation in a standard gel tube with the EBA 270.



hettweb.com/eba270

DISTRIBUTED BY:



Helmer Scientific
 14400 Bergen Boulevard, Noblesville, IN 46060 USA
Toll Free: (800) 743-5637 | **Email:** sales@helmerinc.com
www.helmerinc.com



Discover more solutions online at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
 © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-EBA-270-E Rev A 01-2016.

