

HELMER® SCIENTIFIC Ultra-Low Temperature Freezer Accessories



Racks | Boxes | Dividers | Canisters | Back-up Systems

ULT Freezer Rack Systems


Helmer Scientific offers a wide range of ULT rack systems to safely store and organize your products, samples, and therapies. Designed to last, these storage systems are constructed from high-quality corrosion-resistant stainless steel and are offered in various configurations to categorize and maximize freezer capacity. Options include vial box racks, centrifuge tube racks, microplate racks, and red cell canister racks.

Vial Box Racks



**Vial Box Rack
322098-1**

Rack for 2" vial boxes, side access, 16 box capacity (4 high x 4 deep). Boxes and dividers included.



**Vial Box Rack
322099-1**

Rack for 2" vial boxes, side access, 20 box capacity (5 high x 4 deep). Boxes and dividers included.



**Vial Box Rack
322100-1**

Rack for 2" vial boxes, drawer access, 16 box capacity (4 high x 4 deep). Boxes and dividers included.



**Vial Box Rack
322101-1**

Rack for 2" vial boxes, drawer access, 20 box capacity (5 high x 4 deep). Boxes and dividers included.



**Vial Box Rack
322102-1**

Rack for 3" vial boxes, drawer access, 8 box capacity (2 high x 4 deep). Boxes and dividers included.



**Vial Box Rack
322103-1**

Adjustable rack for 2" or 3" vial boxes, side access. Boxes and dividers included.

| Part No. | Dimensions in/mm wxdxh | Storage | iUF118-GX | iUF126-GX |
|-----------|--|--|----------------------|----------------------|
| 322098-1 | 5.5 x 22.125 x 8.875 140 x 562 x 225 | Boxes per Rack Racks per Compartment Boxes per Freezer | 16 4 320 | 16 6 480 |
| 322099-1* | 5.5 x 22.125 x 11 140 x 562 x 280 | Boxes per Rack Racks per Compartment Boxes per Freezer | 20 4 352* | 20 6 528* |
| 322100-1 | 5.5 x 22.375 x 8.875 140 x 568 x 225 | Boxes per Rack Racks per Compartment Boxes per Freezer | 16 4 320 | 16 6 480 |
| 322101-1* | 5.5 x 22.375 x 11 140 x 568 x 280 | Boxes per Rack Racks per Compartment Boxes per Freezer | 20 4 352 | 20 6 528 |
| 322102-1 | 5.5 x 22.375 x 6.3125 140 x 568 x 161 | Boxes per Rack Racks per Compartment Boxes per Freezer | 8 4 160 | 8 6 240 |
| 322103-1 | 5.5 x 22.375 x 8.875 140 x 569 x 225 | Boxes per Rack Racks per Compartment Boxes per Freezer | 16/8 4 320/160 | 16/8 6 480/240 |

*For use in bottom 2 compartments.

NOTE: Optimized storage in GX ULT includes 4-high racks on top 3 shelves, 5-high racks on bottom 2 shelves.

Centrifuge Tube Racks



**Tube Rack
322104-1**

Rack for storage of 15ml centrifuge tubes, drawer access, 80 tube capacity.




**Tube Rack
322105-1**

Rack for storage of 50ml centrifuge tubes, drawer access, 39 tube capacity.

| Part No. | Dimensions in/mm wxdxh | Storage | iUF118-GX | iUF126-GX |
|----------|--|--|-----------------|-----------------|
| 322104-1 | 5.5 x 22.375 x 4.9375 140 x 564 x 126 | Tubes per Rack Racks per Compartment Tubes per Freezer | 80 4 1600 | 80 6 2400 |
| 322105-1 | 5.5 x 22.375 x 4.9375 140 x 564 x 126 | Tubes per Rack Racks per Compartment Tubes per Freezer | 39 4 780 | 39 6 1170 |

Red Cell Canister Racks



**Vertical Rack
401256-1**

For storage of red cell canisters at ultra-cold temperatures.




**Horizontal Rack
401257-1**

For processing/freezing red cells in canisters at ultra-cold temperatures.

| Part No. | Dimensions in wxdxh | Storage | iUF118-GX | iUF126-GX |
|----------|--------------------------|--|----------------|----------------|
| 401256-1 | 11.2 x 16.3 x 4 | Canisters per Rack Racks per Compartment Canisters per Freezer | 12 2 120 | 12 2 120 |
| 401257-1 | 10 x 15.2 x 4.2 | Canisters per Rack Racks per Freezer | 3 1 | 3 1 |

Microtiter Plate Racks



Microplate Rack 322106-1

Rack for 96 and 384-well microtiter plates, side access, with locking rod.



Microplate Rack 322107-1

Rack for 96 deep-well microtiter plates, side access, with locking rod.



Microplate Rack 322108-1

Rack for 96 and 384-well microtiter plates, drawer access.



Microplate Rack 322109-1

Rack for 96 deep-well microtiter plates, drawer access.

| Part No. | Dimensions in mm wxdxh | Storage | iUF118-GX | iUF126-GX |
|-----------|--|--|---|---|
| 322106-1* | 5.5 x 22.5 x 11 140 x 571 x 280 | Plates per Rack Racks per Compartment Plates per Freezer | 96 w/ lid / 120 w/o lid 4 1920/2400 | 96 w/ lid / 120 w/o lid 6 2880/3600 |
| 322107-1 | 5.5 x 22.5 x 7.5 140 x 571 x 191 | Plates per Rack Racks per Compartment Plates per Freezer | 24 4 480 | 24 6 720 |
| 322108-1 | 5.4375 x 21.75 x 5.3125 138 x 552 x 186 | Plates per Rack Racks per Compartment Plates per Freezer | 54 w/ lid / 72 w/o lid 4 1080/1440 | 54 w/ lid / 72 w/o lid 6 1620/2160 |
| 322109-1 | 5.4375 x 21.75 x 8.0625 138 x 552 x 205 | Plates per Rack Racks per Compartment Plates per Freezer | 24 4 480 | 24 6 480 |

*For use in bottom 2 compartments

ULT Freezer Backup Systems



CO2 Backup System 220628

Protects inventory due to power outage or equipment failure, maintaining temperature down to -60°C using CO2.

- Connects to facility's existing infrastructure for carbon dioxide (CO2) gas delivery
- Injects gas into freezer when temperature threshold has been reached
- System designed to maintain temperature down to -60°C with one 50 lb cylinder of CO2 for a minimum of 10 hours
- Independent display/control module
- Backup battery
- Door safety interlock switch prevents injection while door is open
- Freezer connection hardware
*Does not include gas

LN2 Backup System 220629

Protects inventory due to power outage or equipment failure, maintaining temperature down to -65°C using LN2.

- Connects to facility's existing infrastructure for liquid nitrogen(LN2) gas delivery
- Injects gas into freezer when temperature threshold has been reached
- System designed to maintain temperature down to -65°C with one 65 liter cylinder of LN2 for at least 48 hours
- Independent display/control module
- Backup battery
- Door safety interlock switch prevents injection while door is open
- Freezer connection hardware
*Does not include gas

Vial Boxes | Dividers | Canisters

| Description | Part No. | Quantity |
|--|--|----------|
| Fiberboard Vial Box - 2" (51mm) 5.25 x 5.25" (133 x 133mm) | 322096-1 | 1 |
| Fiberboard Vial Box - 3" (75mm) 5.25 x 5.25" (133 x 133mm) | 322097-1f | 1 |
| Plastic Vial Box - 2" (47mm) 5.12 x 5.12 (130 x 130) Autoclavable molded box w/ transparent lid, 9x9 divider and location ID. | 220620 white 220621 blue 220622 green 220623 pink | 1 |
| Fiberboard Grid Dividers 9x9, 0.5" cell 13mm vials (81) | 220618 | 1 |
| Fiberboard Grid Dividers 10x10, 0.45" cell 12mm vials (100) | 220619 | 1 |
| Red Cell Canister 2000 ml bags 14.8 x 9.2 x 0.8 (376x234x20mm) | 322110-1 | 1 |

Additional Accessories

| Description | Part No. |
|--|--------------------------------------|
| Cryo Gloves Small Medium Large Extra Large | 220624 220625 220626 220627 |
| Chart Paper 4" (102mm) diameter, -100°C to -40°C, pkg of 52 | 800197-1 |
| Remote Alarms Extends alarm range for 1 device/zone or up to 4 devices/zones | 500607-1 RA1 500608-1 RA4 |
| USB Flash Drive | 120640 |
| i.C3° Screen Protectors (pkg of 3) | 450007-1 |